




## Kinetic Stem Cell (KSC) Counting

*Technologies for  
specific quantification and analysis  
of therapeutic tissue stem cells*

# THE CURRENT PROBLEM

Homogeneous



or


Impossible



*How many people think about the dosage of therapeutic tissue stem cell preparations*

# THE CURRENT PROBLEM

Homogeneous  
***Incorrect***



or


Impossible  
***Previously***



*How many people think about the dosage of  
therapeutic tissue stem cell preparations*


# THE ASYMMETREX VISION

*The Reality*



Heterogeneous with difficult to quantify  
low tissue stem cell fraction

# THE ASYMMETREX VISION



*The Asymmetrex Solution:  
Cell Kinetics Distinctions*

# OUTLINE

- I. Tissue stem cell number determines treatment potency
- II. Previous ineffective technologies
- III. Kinetic Stem Cell (KSC) Counting
- IV. KSC Counting Applications

I. Tissue stem cell number determines  
treatment potency

# Tissue Stem Cell Number Is An Important Potency Factor



## **For stem cell therapies (e.g., HSC transplantation medicine)**

Tissue stem cell-specific dosage is a critical quality attribute for:

- Optimizing and reproducing treatment outcomes
- Designing clinical trials
- Normalization for evaluating clinical trial effectiveness

## **For stem cell-engineered gene therapies (e.g., HSC gene editing)**

Monitoring stem cell fraction during engineering

Tissue stem cell-specific dosage is a critical quality attribute

## II. Previous Ineffective Technologies\*

<u>Technology</u>	<u>Effective?</u>	<u>Problems</u>
<b>Flow cytometry</b>	No	Specificity: TSC = CPC
<b>SCID mice</b>	Only for HSCs	\$, mice, 16 wks, unreliable
<b>CFU</b>	No	Specificity: TSC = CPC
<b>Metabolic</b>	No	Specificity: TSC = CPC

\* Except for muscle satellite stem cells (nuclear Pax6, Pax7)


TSC = tissue stem cells

CPC = committed progenitor cells

# III. Kinetic Stem Cell (KSC) Counting


# KSC Counting Principle I

## *In Vivo* TSC Kinetics Continue in Cell Culture




# KSC Counting Principle I

*In Vivo* TSC Kinetics Continue in Cell Culture\*




# KSC Counting Principle II


Therefore, culture total cell output depends on TSC-specific fraction and cell kinetics.



# KSC Counting Begins with Serial Culture



# CPD Curve Culture Dependencies



- Rate
- Maximum
- Time of arrest

## Culture Factors


Input cell number

Split interval

Split fraction

Cell viability


# CPD Curve Cell Kinetics Dependencies



- Rate
- Maximum
- Time of arrest

## Culture Factors

Input cell number  
 Split interval  
 Split fraction  
 Cell viability



## Cell Kinetic Factors

### **TSC Number**


CPC Number  
 TAC Number  
 TSC Viability  
 CPC Viability  
 TAC Viability

### **TSC Symmetric Rate**


TSC CC Times (Asym & Sym)  
 CPC CC Time  
 CPC Division Number

TSC, tissue stem cell; CPC, committed progenitor cell;  
 TAC, terminally-arrested cell; CC, cell cycle;  
 Division Number = number of divisions before producing TACs


# KSC Counting: A computational simulation modeling approach



# The KSC Counting Process




**Serial Cell Count Data**



**TORTOISE Test™ Simulation**

Average of 10 best from 10,000 trials



**KSC Counting Outputs**

- Blue** – Tissue stem cells
- Red** – Committed progenitor cells
- Green** – Terminally-arrested cells



*Note: Analysis for CD34<sup>+</sup> UCB cells*

# KSC Counting

## Rapid Counting Algorithms

For any day of serial culture,  
if you know the culture's PDT,

For any future cell culture,  
if you know the culture's CPD,



...you know the SCF.

...you know the SCF.

**Note:** Data for CD34<sup>+</sup> umbilical cord blood HSCs  
SCF, stem cell-specific fraction

# CD34<sup>+</sup> Selection Validations

[Mean SCF (95CI)]

<b>Source</b>	<b>CD34<sup>-</sup></b>	<b>CD34<sup>+</sup></b>
<b>UCB</b>	0.0001 (0.001-0.001)	0.63 (0.36-0.90)

<b>Source</b>	<b>Unfractionated</b>	<b>CD34<sup>+</sup></b>
<b>UCB-1a</b>	0.03 (0.02-0.05)	
<b>UCB-1b</b>	0.02 (0.01-0.02)	
<b>UCB-2</b>		0.18 (0.14-0.22)
<b>UCB-3</b>		0.72 (0.58-0.85)
<b>MPB-1</b>	0.02 (0.01-0.03)	0.78 (0.61-0.96)
<b>MPB-2</b>	0.15 (0.06-0.24)	
<b>MPB-3</b>		0.83 (0.76-0.91)

**SCF**, stem cell-specific fraction; **95CI**, 95% confidence interval

**UCB**, umbilical cord blood; **MPB**, mobilized peripheral blood

# Human tissue stem cells

## KSC counted to date

- Liver hepatic stem cells
- Lung interstitial stem cells
- Corneal stem cells
- Bone marrow hematopoietic stem cells
- Mobilized peripheral blood hematopoietic stem cells
- Umbilical cord blood hematopoietic stem cells
- Umbilical cord tissue mesenchymal stem cells
- Bone marrow-derived mesenchymal stem cells
- Adipose-derived mesenchymal stem cells
- Oral-derived mesenchymal stem cells  
(bone, gingival, dental pulp)
- Amniotic membrane stem cells

# IV. KSC Counting Applications

## A. Cord Blood Applications

1. Complement CFU for better potency determinations
2. Complement MNC for unit selection for increased yield
3. Complement CFU for more rapid unit stability evaluations

## B. Stem cell therapy and stem cell-engineered gene therapy dosage

## C. Stem cell-based cell and tissue biomanufacturing (*e.g.*, MSCs)

## D. Drug evaluations

1. Tissue stem cell toxicity evaluations to avoid chronic organ failure
2. Regenerative medicine drug evaluations

## E. Environmental chemical analysis for tissue stem cell toxicants (*e.g.*, carcinogens)

## F. Tissue cell research and stem cell biology research

# For more information

Contact:

James L. Sherley, M.D., Ph.D.

President & CEO

[jsherley@asymmetrex.com](mailto:jsherley@asymmetrex.com)

1-617-990-6819

<https://asymmetrex.com/>